

## SmileDiamonds Product Information

### Nothing will Brighten Your Smile Like a Real Diamond!

#### Information on SmileDiamonds:

100% natural diamonds

Cut: Round full cut with 58 facets diamonds

Color: G - I Near Colorless

Clarity: SI1 – I1 or better

SmileDiamonds are high white luster diamonds with brilliant sparkle. The diamonds come mounted ready to be applied.

#### Quick Facts:

- Takes 5 - 10 minutes to be safely and securely applied by a dental professional or a trained reseller
- Surface sealant material is included with the purchase of every SmileDiamond
- Fun & Affordable
- Reversible - Can be removed or replaced
- No Drilling & No Metal
- Bonded to the surface
- No Pain. No damage to the tooth
- Alternative to piercing or tattooing
- 100% natural diamonds
- Cut: Round full cut with 58 facets diamonds
- Color: G - I Near Colorless
- Clarity: SI1 – I1 or better
- Polish: High luster and shine
- No fake CZ's
- It is not a tattoo, it is better than a tattoo.
- Glamorous

\*\*Diamond dimension may vary slightly because all diamonds are hand cut.

<u>Item No.</u>	<u>Base Width</u>	<u>Height</u>	<u>Diamond Diameter</u>
SD-01	3.5mm	.75mm	1.35mm
SD-02	3.5mm	1.15mm	1.65mm
SD-03	3.5mm	1.50mm	1.96mm
SD-04	3.5mm	1.75mm	2.18mm

After the ToothJewelry Sealant has been applied by following the SmileDiamonds Application Instructions around the entire tooth apply one final thin layer over the top of the diamond the finished height is approximately .15mm in addition to the height given above or possible less.

\*\*Once a year or so it is recommended to apply a thin layer of ToothJewelry sealant or Bisco Fortify sealant over the entire SmileDiamond. It's best to do this when you have already setup an appointment with your family dentist.